



wasserLab

AUTWOMATIC PLUS / UV from 3 to 10 l/h

TYPE II WATER

OSMOTIZED WATER



Compact, space-saving device

The Autwomatic Plus range produces and dispenses **Type II Water (ASTM D1193)** and **Osmotized Water** incorporating the latest technology in the production and quality control of purified water, with production versions of 3, 5 and 10 liters per hour with the following quality:

TYPE II WATER

Conductivity	<1 µS/cm
TOC ¹	<30 ppb
Bacterial count ¹	<0.01 cfu/ml
Particles >0.22 µm / ml ¹	<1
Production Flow Rate	3 - 5 - 10 l/h

OSMOTIZED WATER

Conductivity	<98% Rejection from Tap Water
--------------	-------------------------------

1. UV version: These values are typical and may vary depending on the nature and concentration of contaminants in the Tap Water.



Type II Water Applications

- Preparation of culture media
- Preparation of reagents and buffers
- RIA/ELISA
- AA-Flame
- Spectrophotometry

Osmotized Water Applications

- Autoclave feed
- Washer-disinfectant feeding
- Cleaning of laboratory equipment

Versions

MODEL	REFERENCE	FLOW RATE OUTPUT	CONSUMPTION RECOMMENDED
Autwomatic Plus 3	QA03DP	3 l/h	30 l/day
Autwomatic Plus 5	QA05DP	5 l/h	50 l/day
Autwomatic Plus 10	QA10DP	10 l/h	100 l/day
Autwomatic Plus UV 3	QA03DPUV	3 l/h	30 l/day
Autwomatic Plus UV 5	QA05DPUV	5 l/h	50 l/day
Autwomatic Plus UV10	QA10DPUV	10 l/h	100 l/day
TANK OPTIONS ²	REFERENCE		
Pressurized 30 L	70220		
Pressurized 50 L	70230		

2. The system allows the combination of several Tanks.

Water purification process

The equipment integrates various technologies to optimise the water purification process, through the following stages:

Pre-treatment

The equipment is designed with a pre-treatment system to guarantee the protection of the osmosis membrane, eliminating particles smaller than 1 micron, which contributes to the reduction of mineral incrustation, organic matter and the elimination of chlorine.

The depth filter presents a high filtration capacity, with optimised re-

tion of the colloids present in the water.

The granular bacteriostatic activated carbon is effective in removing free chlorine and minimising bacterial growth.

The polyphosphate-based anti-limescale agent protects against scaling, preventing the precipitation of calcium and magnesium salts inside the equipment, without releasing ions.

The 10 l/h equipment is designed with an external module in accordance with its higher production capacity.



Committed to the environment: EFFICIENT SOLUTIONS THAT SAVE WATER AND ENERGY

Reverse Osmosis

The high performance reverse osmosis system removes up to 99.95% of the organic matter present in the water and up to 98% of the Total Dissolved Solids (TDS). In addition, the equipment has an automatic washing system, designed to prolong its useful life.

Storage¹

The water produced in the reverse osmosis phase is accumulated in a hermetic, pressurized tank of 30 or 50 liters, being protected from air and environmental pollution.

The system allows to manually dispense Osmotized Water from a stopcock placed in the tank outlet, or to automatically feed washer-disinfector or autoclaves. Connections to these types of systems are not included in the installation kit of this equipment, so it is recommended to consult the factory.

1. The Autwomatic Plus equipment can be configured to store Type II Water in atmospheric Tanks.



Purification Phase with Ion Exchange Resins

The Osmotized Water passes through a mixed bed cationic/anionic ion exchange resin, retaining the few dissolved salts in the water, providing purified Type II Water Analytical Grade, according to ASTM (American Society for Testing and Materials) specification, with a conductivity of less than 1 µS/cm.


The configuration of the system by means of a Pressurized Tank of Osmotized Water plus a subsequent ion exchange phase **provides, unlike other systems, a freshly produced Type II Water.** In other words, the user does not store Type II Water but obtains the water directly from the equipment, with a sufficient flow rate to meet their needs.

UV lamp and 0.22 µm final filter (UV version)

To guarantee the microbiological control of purified water, the equipment is equipped with an Ultraviolet lamp that performs bacteriostatic and germicidal functions. This lamp emits a wavelength of 254 nm, which makes it effective in eliminating micro-organisms without compromising water quality.

To meet even more stringent microbiological requirements (<0.01 cfu/ml), the equipment incorporates a 0.22 µm in-line final filter. This additional filter ensures the retention of possible microorganisms that may be present in the final water, guaranteeing that the purified water meets high microbiological quality standards.

Hydraulic system diagram Autwomatic Plus / UV



Tap water



Pressure regulator



Module of pre-treatment



Probe conductivity



Pump



Reverse Osmosis Membrane



Pressurized Tank



Module of deionization



UV lamp






Resistivity and temperature probe



Final filter 0.22 µm

Operation and monitoring



1. Dispensing

The equipment is designed to operate automatically, ensuring that the tank is always kept full, thanks to its automatic stop system. In addition, its design ensures simple and accessible operation for the user.

It incorporates a 4.3 inch touch screen, which facilitates the dispensing of purified water in various ways, adapting to the user's needs.

The available options are:

- Continuous dispensing
- Dispensing by volume¹
- Dispensing by time

1. Accuracy not suitable for flush volumes.

2. Monitoring

This comprehensive monitoring system allows detailed tracking of key aspects of the process, ensuring that the equipment is operating within ideal parameters and guaranteeing the quality of the water produced.

- Conductivity measurements (at 25°C):
 - Equipment feed water.
 - Permeate water from the reverse osmosis module.
 - Type II Water produced.
- Percentage of performance of the reverse osmosis module.
- Control of operating parameters:
- Working hours of the different components of the system.
- Liters produced during the purification process.

3. Customization and Security

The system offers Customization options, allowing to adjust the type of dispensing and conductivity settings according to the user's needs. In addition, to ensure security and access control, the equipment has a custom password that allows the user to access specific menus and functions, ensuring that only authorised persons can make adjustments or modifications to the system.

hh:mm	dd.mm.yyyy	Resetting counters		Pumps	
	Time (h)	Working	Reset	P61	350
Pretreatment	500	500	Stop	P21	60
UV/Photooxid	500	500		P24	0
Final filter II	350	350		P46	0
Type II resin	350	350	Total hours	0	
Osmosis 1	90		Check equip.	0	
			L. Tot II	0	L. Tot I 0
			Pass		Exit

hh:mm	dd.mm.yyyy	wasserLAB	
		19,0° C	Type II water
			Menu

4. Automatism

The system is equipped with a microprocessor that manages the automatic start and stop of the equipment, depending on the volume of water accumulated in the tank. It also has various automated systems to guarantee optimum performance and prolong the useful life of the equipment, such as:

- Automatic shutdown in the event of a water cut-off, to prevent operation without supply.
- Cleaning of the osmosis membrane by flushing, which helps to maintain the efficiency of the filtration process.
- Cleaning of the osmosis membrane with osmotized water, in order to prolong its useful life and ensure maximum performance.

These automatism contribute to an efficient and low maintenance operation, ensuring the constant quality of the purified water.

5. Data Output

The equipment is designed to allow the extraction of operating data to an external memory (USB). The report generated includes detailed records on the quality and quantity of water dispensed, as well as maintenance alerts and changes made to consumables, providing a useful tool for monitoring and controlling system performance.



6. Mobile phone notification (optional)

The system can send alarm notifications directly to mobile phones, allowing real-time alerts on various problems or irregularities in the operation of the equipment.



**Easy and efficient
maintenance:
A SIMPLE AND
FAST PROCESS**

Preventive maintenance, sanitization and calibration

Ease of maintenance and control of the system

The system has been carefully designed for ease of maintenance, allowing the user to perform tasks easily and efficiently. The replacement of consumables is carried out quickly, thanks to a quick-connect system with anti-drip technology built into the cartridges.

The service life of the consumables depends on several factors, such as the quality of the incoming water, including its turbidity, hardness and conductivity, as well as the amount of water dispensed over time.

The integrated software is configured to perform scheduled self-checks, ensuring constant and effective monitoring of the system's operation. This control ensures continuous monitoring of the equipment components and the values related to the quality of the water produced.

In addition, the system issues warnings to notify the user about the need to change consumables, water cuts or possible malfunctions of the measuring probes, allowing an early intervention in case of incidents.

System Sanitization

The system is designed to facilitate the sanitization of the equipment through a semi-automatic process, which ensures a thorough and effective cleaning of all its components. During sanitization, the equipment performs a series of automated steps that include the circulation of disinfectant solutions through critical parts of the system, such as membranes and filters. User intervention is limited to initiating and monitoring the process, following clear instructions provided by the system. This sanitization process is designed to remove microorganisms, bacteria and other contaminants that may have accumulated in the equipment, ensuring that the system continues to operate at maximum efficiency and that the water produced always maintains the highest quality standards. The function also helps to extend the life of the equipment by preventing the build-up of impurities that can affect its performance.

Additional system functionalities

Optional Remote Dispensers

Digitally controlled remote dispensers designed to allow additional Type II water outlets at a distance from the main equipment, optimising space and improving operational efficiency.

Pedestal support option (Ref. W-DIS100-A03) Wall option (Ref. W-DIS102-A03)



Equipment can be integrated into furniture

The equipment is designed to be fully integrated into laboratory furniture, optimising the available space and leaving the laboratory work bench free for other tasks. Its minimalist design adapts perfectly to laboratory work environments, offering an aesthetic and functional solution that maximises efficiency without compromising system performance. We work with leading furniture brands.



Manual Dispensers

Manual Dispenser (Ref. W-DIS006)

Additional mechanically operated osmotized water outlet, especially convenient for filling carafes and dispensing a few meters away from the main equipment.

Flexibility to offer solutions that ADAPTED TO EACH LABORATORY

Accommodating the needs of the space available

Wall bracket (Ref. 10261)

Base designed to allow safe and stable installation of the equipment directly on the wall. Its robust structure guarantees a firm mounting, optimising the use of space and ensuring that the equipment is well fixed and accessible. Ideal for environments where space in the work area needs to be freed up.

Compact Module (Ref. 10092)

A functional and compact design cabinet, it offers a solution for housing the equipment and its components in an orderly and efficient manner. Perfect for environments where equipment needs to be kept protected and in place, while ensuring accessibility and ease of use.



Wall bracket

Compact Module

Design and Installation of Distribution Loops

We design and install distribution loops, interconnected systems that guarantee an efficient distribution of purified water between different points, adapted to the specific needs of each project.

IQOQ Qualifiable Equipment for the Pharmaceutical Sector

The equipment is designed to be qualified in the processes of IQOQ (**I**nstallation and **O**perational **Q**ualification) required in the pharmaceutical sector. It complies with industry specific regulatory standards, ensuring its suitability for use in regulated environments, where traceability, quality and process validation are critical to ensure compliance with current regulations.

Declaration of Product Use: WEEE Directive

In accordance with European Union legislation, this product will be considered **Waste Electrical and Electronic Equipment (WEEE)** once it reaches the end of its useful life.

For detailed information on the recycling and proper disposal of this product, please contact our website.

Quality Assurance to Facilitate GLP and cGMP Compliance

The system has been designed and manufactured to facilitate its integration into regulated working environments such as GLP (Good Laboratory Practices) and cGMP (current Good Manufacturing Practices). Some of its outstanding features include:

- Manufactured under the standards ISO 9001:2015 and ISO 14001, ensuring that the product meets the highest standards of quality and environmental management.
- **CE marking:** The equipment has passed rigorous safety and electromagnetic compatibility tests (emission and immunity), carried out by an external accredited centre, which certifies its compliance with European safety and performance standards.
- **Calibration certificate:** The equipment is delivered factory calibrated, guaranteeing its accuracy from the first moment of use. It also allows the adjustment and recalibration of the conductivity meter by means of a certified standard, traceable to the national standards of the Deutscher Kalibrierdienst (DKD) of Germany, ensuring the reliability and accuracy of the measurements over time.



Installation Requirements

- Alternating current socket 110 / 120 / 230 V - 50 - 60 Hz. with earth connection at a maximum of 2 meters from the equipment.
- Tap water connection (maximum 3 meters).
- Connection: 3/8" male gas thread.
- Drainage (maximum 3 meters).
- Tap water quality:
 - Conductivity: <math><2.000 \mu\text{S} / \text{cm}</math>
 - pH: 4 - 10
 - Hardness: <math><300 \text{ ppm CaCO}_3</math>
 - Turbidity: <math><1 \text{ NTU}</math>
 - CO_2 : <math><30 \text{ ppm}</math>
 - Silica: <math><30 \text{ ppm}</math>
 - TOC: <math><50 \text{ ppb}</math>
 - Free chlorine: <math><1.5 \text{ ppm}</math>
 - SDI: <math><7</math>
 - Temperature: 5 - 35 °C
- Pressure: 2 - 6 bar.
- Installation space for the equipment and its elements, guaranteeing an accessible work area for handling.

Specifications:

Dimensions:

- Equipment Autwomatic Plus / UV: 60 x 36 x 49 cm (height / width / depth).
- External pre-filtration (10 l/h model): 55 x 24 x 16 cm (height / width / depth).
- 30 liters tank: 60 cm height x 40 cm diameter.
- 50 liters tank: 80 cm height x 40 cm diameter.
- Compact module: 96 x 46 x 60 (height / width / depth).

Weight: 35 Kg.

Power consumption: 0.6 A (230 VAC) - 1.2 A (110 VAC).


Power: 136 VA (230 VAC) - 136 VA (110 VAC).

Noise level: <math><50 \text{ dB}</math>.


Communication port: USB.

Language Software: English, Spanish, French, Portuguese and Italian.

Equipment



Tanks





Wasserlab

Water Purification Systems

Wasserlab®

We are manufacturers of **water purification equipment** with an extensive track record in the installation of solutions in **multiple sectors**.

We offer **personalised advice** in the selection of equipment and we provide **comprehensive technical support** to guarantee optimum operation.

Pol. Ind. Comarca II Calle E, Nº 3. 31191 Barbatáin - Navarra (Spain)
T. +34 948 186 141 - info@wasserlab.com
www.wasserlab.com

